SOMMAIRE

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p 6

GENERAL INSTRUCTIONS
EINSATZ DES GERÄTS
UTILIZZO DELL'APPARECCHIO
UTILIZACIÓN DEL APARATO
ИСПОЛЬЗОВАНИЕ АППАРАТА
ANVÄNDA TRÄNINGSPAPPARATEN
UTILIZAÇÃO DO APARELHO
GEBRUIK VAN HET APPARAAT

II- LA FONCTION ELECTRODES

p 8

ELECTRODES
DIE ELEKTRODENFUNKTION
LA FUNZIONE ELETTRODI
LA FUNCIÓN ELECTRODOS
ФУНКЦИЯ ЭЛЕКТРОДОВ
ANVÄNDNING MED ELEKTRODER
A FUNÇÃO ELÉTODOS
DE ELEKTRODENFUNCTIE

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p 9

BELT
DIE GÜRTELFUNKTION
LA FUNZIONE CINTURA
LA FUNCIÓN CINTURÓN
ФУНКЦИЯ ПОЯСА
ANVÄNDNING MED BÄLTE
A FUNÇÃO CINTA
DE GORDELFUNCTIE

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p 9


ADVICE ON USE
EINSATZSRATSCHLÄGE
CONSIGLI PER L'USO
CONSEJOS DE UTILIZACIÓN
РЕКОМЕНДАЦИИ ПО ИСПОЛЬЗОВАНИЮ
ANVÄNDNINGSRÅD
CONSELHOS DE UTILIZAÇÃO
GEBRUIKSADVIEZEN

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p 11


TECHNICAL CHARACTERISTICS
TECHNISCHE DATEN
CARATTERISTICHE TECNICHE
CARACTERÍSTICAS TÉCNICAS
ТЕХНИЧЕСКИЕ ХАРАКТЕРИСТИКИ
TEKNISK SPECIFIKATION
CARACTERÍSTICAS TÉCNICAS
TECHNISCHE EIGENSCHAPPEN

SPORT-ELEC BODY BEAUTIFUL



IV

FONCTION ELECTRODES

V

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VI

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INTRODUCTION TO ELECTRO-MUSCULATION

What does it do ?

Electro-musculatation allows muscles to be developed and strengthened and to improve in condition, tonus, preparation for effort and recovery after an effort.

The aim is to maintain the beauty of the body, make sport easier and to keep muscles in good health!

Depending on our activity, certain muscles can be more or less developed, while others, not used so much, tend to shrink.

Our organism consists of over 600 muscles. Bringing all the muscular fibres of a muscle into use requires intense and prolonged physical efforts. In electro-stimulation, this work is performed without effort. The fitness and harmony of our body depends on our muscles. Electro-musculatation using **SPORT-ELEC** allows you to work your muscles whatever your physiology.

The quality of **SPORT-ELEC** currents :

Currents are defined according to several parameters :

Intensity: measured in milliamperes (ma) and determines the extent of muscular contraction. The intensity required for complete contraction depends on the individual.

Frequency : this is measured in Hertz (Hz) and represents the number of pulses per second. It determines the type of muscular fibres worked. Endurance work will require a frequency of around 20 Hertz, while explosive muscular work will require a frequency of around 100 Hertz. .

Pulse width : this is measured in microseconds (μ s) and makes it possible to modulate the strength of stimulation. Muscles do not all have the same receptivity, and some can be difficult to arouse due to their lack of use. These therefore require a higher pulse width.

"VARIABLE FREQUENCY" :

SPORT-ELEC machines are equipped with powerful microcontrollers, able to provide programmes made up of varying frequencies. Alternate, scanned (or staggered) and set frequencies are also possible. The more frequencies vary, the more comfortable sessions are.

How does it work ?

Electro-musculatation is the action which generates a contraction of the muscular fibres.

Natural contraction of the motor muscles is caused by electrical pulses coming from the central nervous system.

Electro-musculatation causes a muscular contraction obtained through electrical pulses coming from outside (power generators) through electrodes placed in contact with the skin, on the muscle to be contracted.

Who is it aiming at ?

All those wishing to maintain the good condition of their body and the health of their muscles.

SPORT-ELEC is designed for all types of muscle form.

What are the limitations of its use?

The only contra-indications concern people with pacemakers and pregnant women (see "safety instructions/ contra-indications") and application over electronic implants.

Recommendations

It is essential, before using any electro-stimulation appliance, to read and understand its instructions for use.

Results can be visible after several sessions and may vary depending on the musculature of individual users.

If you are suffering from a chronic illness, you are recommended to consult a doctor.

GENERAL

Body Beautiful, thanks to its connections, can link the multi-function belt with CEV₁ (photo V-a page 4) for the abdomen, thighs, buttocks and the back **and** with CEV₂ (photo IV-b page 3) to the adhesive electrodes to work muscles or areas of muscle precisely on all parts of the body.

Body Beautiful has 2 independent generators and 4 programmes including 35 working phases, to be used with the belt or with the adhesive electrodes.

Presentation of Body Beautiful (photo I, page 2)

- | | |
|--|---|
| a - Belt | d - Adhesive electrodes (90 x 50 mm) |
| b - Machine | e - Machine/belt connecting strap (CEV1) |
| c - Humidifiable long-duration electrodes (100 x 50 mm) | f - Machine/electrodes connector (CEV2) |

I. GENERAL INSTRUCTIONS

Whether you use the "belt function" or the "electrodes function", the start procedure for your Body Beautiful is identical.

Presentation of the machine (photo II, page 2)

- | | |
|---|--|
| a - Diodes indicating the current programme | e - Button for grey power setting (Channel B) |
| b - Images indicating programmes | f - On/Off |
| c - Programme selection button | g - Connection of connectors CEV1 and CEV2 |
| d - Button for black power setting (Channel A) | |

- 1°** Check that the switch for the **Body Beautiful** is on "0" and that the power-on buttons are at minimum (pushed to the left). Open the "battery" compartment and insert 2 x 9 V batteries into the compartment and then close the case. Comply with + and - indications.

IMPORTANT

(see diagram III page 2)

INSERTION OF BATTERIES : connect the tags to the batteries and insert the back of the battery by pressing with the thumb. Follow these instructions in the order given to avoid damaging the tags.

- Gripper clip
- Battery compartment
- 9 volt batteries (not provided)
- Detachable cover

2° START OF SESSION / PROGRAMME SELECTION

a) Once you have selected the electrode function (see page 8), or the belt function (see page 9), you can start the session. Make sure that the two buttons for the Black (channel A) and Grey (channel B) channels are at minimum (safe position for starting up the appliance). Move switch to "I".

b) CHOOSE YOUR PROGRAMME :

Simply press on the "programme selection" button under the screen. The programmes will appear. Programmes are shown by schemas. When a programme is selected, the corresponding diode above the schema lights up.

INVIGORATING MASSAGE

20 mins prior to sport to prepare and warm up your muscles, or to relax after making an effort.

Staggered frequencies vary from 1 to 13 Hz to provide a genuine massage effect.

Total Duration : 20 min
Repeats : 47

SILHOUETTE SCULPT

40 mins of body building to firm up, but also to increase muscular volume after a period of inactivity. This programme is ideal for keeping the body firm and for re-shaping your silhouette.

Frequencies vary from 3 to 60 Hz alternated with frequencies at 10 Hz for active recovery to build up muscle resistance.

Total Duration : 40 min
Repeats : 34

MUSCLE TONIC

20 mins of more strenuous weight lifting to work on target areas.

Frequencies vary from 70 to 100 Hz alternated with frequencies at 3 Hz for active recovery to build up tonicity and muscular strength.


Total Duration : 20 min
Repeats : 27

RELAXING MASSAGE (TENS)

20 mins to relieve muscular fatigue and allow you to recovery after physical effort.

Staggered frequencies vary from 40 to 100 Hz to provide the state of wellbeing of a relaxing massage.

Total Duration : 20 min
Repeats : 63

3° Intensity setting

Gently turn the black button (channel A) from left to right, you will begin to feel contractions (from the level "2"), increase the power if necessary. Contractions should be firm but never painful.

Continue in the same way for the grey power button (channel B).

If you feel an unpleasant sensation around the electrodes: bring the power-on buttons back to minimum, move the electrode(s) and gradually increase the power. Once the adjustments have been optimised in this way, leave your appliance to work for the full duration of the programme.

- 4°** Once the programme has ended, return the power-on buttons to minimum and move the switch back to the "0" position.

— II. ELECTRODES FUNCTION —

This function of **Body Beautiful** allows you to use the 4 adhesive electrodes by connecting them to the machine.

1° Take the 4 electrodes and carefully remove the holder which protects each of them, moisten them one by one, spreading a few drops of water on the adhesive gel with the tips of your fingers and wait for 2 minutes. (see **diagram IV-a page 3**)

2° Position the electrodes on your body as shown in diagrams A, B, C, D, E, F, G, H **page 3** :

Examples of use : A = channel A (black button) ; B = channel B (grey button)

BE CAREFUL: Please place the electrodes horizontally inside the belt (see diagram c page 4)

A – Pectoral muscles	4 – Back	9 – Upper trapezius	15 – Outer gemelii
1 – Major pectoral	C – Abdominal muscles	10 – Middle trapezius	16 – Inner gemelii
a – Do not place electrodes in the heart region	5 – Major oblique	F – Back and buttocks	H – Thighs
b – Heart	6 – Major longitudinal	11 – Main back muscles	17 – Major adductors
B – Shoulders	D – Biceps	12 – Main gluteal muscles	18 – Vastus internus
2 – Middle deltoid	7 – Biceps - longitudinal	G – Legs	19 – Rear longitudinal
3 – Front	8 – Biceps - transversal	13 – Crural biceps	20 – Vastus externus
	E – Trapezi	14 – Half tendinous	

3° Take the CEV2 connector (photo **IV-b** page 3) and plug it into the machine (photo **IV-c** page 3). Then take the 4 electrodes and connect these plugs to the plugs of the connecting wires (see photo **IV-d** page 3).

4° The box is connected to the electrodes via the CEV2 connector and by two groups of 2 80cm-long wires ending in plugs of different colours (2 black and 2 grey).
You may attach your **SPORT-ELEC** to the inside of an article of clothing thanks to the clip.

5° Turn your machine on following the procedure under **“GENERAL INSTRUCTIONS”** as described on page 13 of this brochure.

6° Once the programme has ended, return the power-on buttons to minimum and move the switch back to the **“0” position**. Gently disconnect and remove the electrodes, moisten them slightly with a few drops of water and put them carefully back on their holder. Electrodes are stored in their original packaging in a cool place (vegetable compartment of a fridge).

— III. BELT FUNCTION —

1° PRESENTATION OF THE BELT :

Front side : **photo V-a page 4**

1- Belt

2- CEV1 strap

3- Belt setting zone

Electrode side : **photo V-b page 4**

- 8 female plugs to attach the 4 humidifiable electrodes

The belt is made up of a nylon section and a thermoformed section. Total belt size: 98 cm.

The belt is adjustable from 34 to 56.

2° MACHINE CONNECTION TO THE BELT

To use your **SPORT-ELEC** belt :

1 - Place the belt on a flat surface, the female plugs facing you.

2 - Check that your **SPORT-ELEC** machine is on OFF.

Connect the machine to the CEV1 strap (**photo V-d page 4**)

3 - Pass the yellow side of the electrodes (**photo I-c page 2**) under warm tap water in order to wet them.

Drain them slightly to avoid the water dripping excessively.

4 - Attach the electrodes to the belt using the plugs (voir **photo V-b page 4**).

5 - Place the belt so that it lies centrally on the abdomen (**schéma V-c page 4**).

Do not hesitate to place the belt slightly to the left or right (check that your **SPORT-ELEC** is OFF).

The setting zone in the centre of the belt allows you to adjust the electrodes depending on the muscular zone. You can tighten or loosen the belt using the drawstring in the middle of it.

3° BELT START

WARNING : Do not turn your **SPORT-ELEC** on before attaching the belt.

Once your machine is connected, turn it on and following the **GENERAL INSTRUCTIONS** on page 13.

4° BELT STOP

Once the programme has ended, turn the power buttons to minimum and turn the machine OFF. Detach the machine from the strap.

— IV. ADVICE ON USE —

CONTRA-INDICATIONS :

Do not use **SPORT-ELEC Body Beautiful** :

-On the cardiogenic area

-If you have been fitted with a pacemaker

-If you are pregnant

SAFETY PRECAUTIONS :

* Do not use body electrodes on the face

* Do not leave within the reach of children

* Ask for advice from your doctor :

- In case of electronic or metal implants in general (pins, etc..)

- In case of dermatitis (wounds, etc..)

* The appliance is not recommended for therapeutic use (seek advice from your doctor).

WARNING

Do not use your **SPORT-ELEC Body Beautiful** in a damp atmosphere (bathroom, sauna) or close to an electrical appliance connected to the mains and the ground, or even piping connected to the ground.

A patient's simultaneous connection to high-frequency surgical equipment may cause burns at the contact points of the stimulator's electrodes and the stimulator may possibly be damaged. Operation in the immediate vicinity of a short-wave therapy apparatus may cause instabilities in the power output from the stimulator.

IDEAL USE BETWEEN SIX WEEKS AND THREE MONTHS AFTER GIVING BIRTH

to help you find your muscular tonus once again

USER TIPS

- It is possible to restart the cycle of one of the programmes on another muscular zone without however exceeding the 40 min of use on the same part of the body.

If you wish to stop a programme in progress, all you have to do is put the 2 buttons back to "o" and repeat the procedure for starting up the appliance. Then select a new programme and repeat the procedure for starting up the appliance.

- For optimal use, take inspiration from the diagrams reproduced in this brochure.

- The case and wires of SPORT ELEC may be cleaned with a blotter moistened with a disinfectant solution (such as Mercryl) diluted x 10 with water]

- Parts in contact, such as electrodes, may be rehydrated before and after use with a disinfecting agent (e.g. Hèxomédine)

- The electrodes must be repositioned on their holder and stored in a cool place between +5°C and +10°C.

- The electrodes have a life of approximately 40 applications, don't hesitate to order them from us directly at the company's head office.

SPORT-ELEC S.A.,

B.P. 35 - 27520 BOURGTHEROULDE

TEL. : + 33 (0)2.32.96.50.50

or contact your local supplier

SAFETY CHECK

Every 24 months, the following points must be checked:

1- Check whether there has been superficial deterioration of the material

2- Check that the appliance is functioning correctly:

* legibility of writing on the case

* switch, contactor

3- Check the reliability of accessories:

* wire from electrodes in good condition
* electrodes in good condition (no cracks, no defects in material)

4- Measure electrical values:

* for a charge of 500 ohms, the peak voltage at outlet must not exceed 100 V.

ENVIRONMENTAL CONDITIONS FOR TRANSPORT AND STORAGE :

- Ambient temperature between +10°C and + 40°C

- Relative humidity between 30% and 75%

- Pressure 700 hPa

BELT MAINTENANCE AND CLEANING ADVICE

- WARNING: Do not put in washing machine or immerse in water.

- Do not dry clean

- Do not use bleach

- Clean the belt with a humid sponge with detergent. Wipe with a cloth.

- Clean the box with a cloth with disinfectant (such as Mercryl) diluted with water to 1/10e.

- Clean the electrodes with a mixture of water/hypoallergenic detergent, rinse well.

- Store the adhesive electrodes in a cool location in a plastic bag (the vegetable rack in the fridge for example).

- Please dry rivets of the belt and electrodes after use.

Recommendations :

The machine, electrodes and belt may be disposed of with environmentally friendly household waste.

The used batteries should be disposed of in existing collecting systems or special containers.

V. TECHNICAL CHARACTERISTICS OF BODY BEAUTIFUL

Current type	Two phase
Number of programmes	4
Total number of phases	35
Number of phases per programme	6 - 60
Number of electrodes	4 adhesive (90 x 50 mm) 4 humidifiable long-duration (100 x 50 mm)
Number of adjustable independent generators	2
Maximum intensity for 500 Ω / channel	 86 mA / 500 Ω / channel
Frequency range	From 1 to 100 Hz
Pulse range	From 160 to 320 µs
Power supply	2 9 V batteries (not supplied)
Type of batteries	6LR61 - 6F22
Battery life cycle	1 month for 20 minutes use per day
Dimensions:	
- machine	75 x 145 x 40 mm
- belt (adjustable from 34 to 56)	970 mm long
Total weight	558 g
Automatic stop at end of programme	Yes
Start safety	100%
BF classification	
Warning : read instructions carefully	
Standards	 (89/336/CEE) CEM
Warranty	2 years
Made in France	Yes

Body Beautiful is a registered trademark. The box is in polished ABS.

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A envoyer à **SPORT-ELEC S.A.**
Fabriqu  en France par / Made in France by **SPORT-ELEC S.A.**